



Bottles, glass, powder, extra wide screw neck

02

Clear, with black plastic screw cap and wood pulp faced steran liner. (PVDC).
For spare caps, see page 181.

Catalogue No	Capacity, mL	Thread style	Height with cap, mm	O.D., mm	Pack qty
BTF-625-010X	15	R3/38	37	41	40
BTF-625-030Y	30	R3/38	48	41	20
BTF-625-050C	60	R3/51	55	52	20
BTF-625-070T	120	R3/58	65	62	20
BTF-625-090N	250	R3/83	75	86	10
BTF-625-110K	500	R3/100	86	110	10



Bottles, glass, powder, extra wide screw neck

02

Clear, with fitted white polypropylene cap without liner.
Caps are ideal for use with marker pen. For spare caps, see page 181.

Catalogue No	Capacity, mL	Thread style	Height with cap, mm	O.D., mm	Pack qty
BTF-630-030Y	30	R3/38	48	41	20
BTF-630-050L	60	R3/51	55	52	20
BTF-630-070T	120	R3/58	65	62	20
BTF-630-090N	250	R3/83	75	86	10



Fisher Scientific Chemicals Catalogue 2009/10

Analytical products, services and support from the market leader
in the supply of laboratory pack analytical chemicals.