

Bottle carrier, 4-in-1, Thermo Scientific Nalgene®, EZ Tote™



54

- Versatile carrier holds 500mL, 1L or 4L bottles
- Easy to assemble and stores flat
- Assembly instructions are printed on each tote
- Good chemical resistance
- Assembled weight 90g



Catalogue No	Alt. No	Dimensions [w x d x h], mm	Pack qty
BTR-846-010J	6565-0001	200 x 330 x 343	4

Bottle carriers, polyethylene, Winchester

02



Accommodates one or two 2.5L Winchester shape bottles. Liquid tight, with carrying handle pivoted for easy pouring.

Catalogue No	Alt. No	Type	Colour
BTS-850-S	BGS 012	Double	Red
BTS-852-G	BGS 014	Double	White

Bottle carrier, polythene

02



Holds 6 x 500mL bottles.

Catalogue No	Alt. No	Diameter, mm
BTS-861-020W	BGS 024	90 (compartment)

Bottle carrier, polypropylene, wash

03



Blue carrier features a secure handle and wells for carrying up to 6 x wash bottles

Autoclavable at 121°C and chemical resistant. Holds either 250mL or 500mL bottles. Hole diameter 73.152mm.

Catalogue No	Alt. No	Dimensions [w x d x h], mm
BTS-920-030W	F11640-0001	286 x 194 x 225