

Test tube racks, Thermo Scientific Nalgene[®], ResMer[™], Unwire[™]

55




a

Made of tough chemically resistant ResMer[™], they hold a full load of tubes without twisting or bending. Moulded in one piece with a smooth finish and fewer places for contamination to lodge. Excellent thermal properties and can be repeatedly autoclaved up to 121°C. Does not get brittle or stick to fingers when at sub-zero temperatures. Large end plates permit clear labelling. Letters and numbers are moulded into top tiers for easy, accurate identification of individual tubes. Racks of the same size stack when empty.

Catalogue No	Alt. No	Dimensions [w x d x h], mm	Tube diameter [maximum], mm	Colour	Array
STK-770-010A	5970-0013	102 x 200 x 57	13	White	6 x 12
STK-770-012T	5970-0113	102 x 200 x 57	13	Orange	6 x 12
STK-770-014P	5970-0213	102 x 200 x 57	13	Yellow	6 x 12
STK-770-016L	5970-0313	102 x 200 x 57	13	Blue	6 x 12
STK-770-018H	5970-0413	102 x 200 x 57	13	Green	6 x 12
STK-770-020U	5970-0513	102 x 200 x 57	13	Red	6 x 12
STK-770-030R	5970-0016	127 x 248 x 70	16	White	6 x 12
STK-770-032N	5970-0116	127 x 248 x 70	16	Orange	6 x 12
STK-770-034J	5970-0216	248 x 127 x 70	16	Yellow	6 x 12
STK-770-036F	5970-0316	248 x 127 x 70	16	Blue	6 x 12
STK-770-038B	5970-0416	248 x 127 x 70	16	Green	6 x 12
STK-770-040X	5970-0516	248 x 127 x 70	16	Red	6 x 12
STK-770-050L	5970-0020	250 x 102 x 83	20	White	4 x 10
STK-770-052H	5970-0120	250 x 102 x 83	20	Orange	4 x 10
STK-770-054D	5970-0220	250 x 102 x 83	20	Yellow	4 x 10
STK-770-056W	5970-0320	250 x 102 x 83	20	Blue	4 x 10
STK-770-058S	5970-0420	250 x 102 x 83	20	Green	4 x 10
STK-770-060Y	5970-0520	250 x 102 x 83	20	Red	4 x 10
STK-770-070F	5970-0025	300 x 121 x 92	25	White	4 x 10
STK-770-072B	5970-0125	300 x 121 x 92	25	Orange	4 x 10
STK-770-074U	5970-0225	121 x 300 x 92	25	Yellow	4 x 10
STK-770-076Q	5970-0325	121 x 300 x 92	25	Blue	4 x 10
STK-770-078M	5970-0425	121 x 300 x 92	25	Green	4 x 10
STK-770-080C	5970-0525	121 x 300 x 92	25	Red	4 x 10
STK-770-090W	5970-0030	108 x 283 x 83	30	White	3 x 8
STK-770-092S	5970-0130	108 x 283 x 83	30	Orange	3 x 8
STK-770-094X	5970-0230	108 x 283 x 83	30	Yellow	3 x 8
STK-770-096K	5970-0330	108 x 283 x 83	30	Blue	3 x 8
STK-770-098G	5970-0430	108 x 283 x 83	30	Green	3 x 8
STK-770-110T	5970-0530	108 x 283 x 83	30	Red	3 x 8

Test tube racks, half rack, Thermo Scientific Nalgene[®], ResMer[™], Unwire[™]

55



a

All the features of the full rack (STK-770 series above), but for fewer tubes.

Catalogue No	Alt. No	Dimensions [w x d x h], mm	Tube diameter [maximum], mm	Colour	Array
STK-780-030K	5972-0016	126 x 126 x 68	16	White	6 x 6
STK-780-036C	5972-0316	126 x 126 x 68	16	Blue	6 x 6
STK-780-038R	5972-0416	126 x 126 x 68	16	Green	6 x 6
STK-780-040H	5972-0516	126 x 126 x 68	16	Red	6 x 6
STK-780-050E	5972-0020	128 x 103 x 83	20	White	4 x 5
STK-780-058L	5972-0420	128 x 103 x 83	20	Green	4 x 5
STK-780-070V	5972-0025	122 x 122 x 75	25	White	4 x 4
STK-780-076J	5972-0325	122 x 122 x 75	25	Blue	4 x 4
STK-780-090P	5972-0030	109 x 109 x 84	30	White	3 x 3
STK-780-096D	5972-0330	109 x 109 x 84	30	Blue	3 x 3
STK-780-098W	5972-0430	109 x 109 x 84	30	Green	3 x 3
STK-780-100P	5972-0530	109 x 109 x 84	30	Red	3 x 3