

Test tube racks, nylon coated

03

a


Catalogue No	No. of holes	Dimensions [w x d x h], mm
STK-750-010X	12	200 x 78 x 60

Nylon coated steel. Autoclavable to 130°C and with good resistance to chemicals and abrasions

Catalogue No	No. of holes	Tube diameter [maximum], mm	Tube length [maximum], mm
STK-640-030A	36	10 to 14	50
STK-640-050R	12	13 to 18	100
STK-640-070L	24	13 to 18	100
STK-640-090F	18	13 to 21	100
STK-640-110C	36	13 to 21	100
STK-640-130T	12	16 to 25	100

Nylon coated steel. Autoclavable to 130°C

Catalogue No	Hole diameter, mm	Dimensions [w x d x h], mm
STK-755-010H	18 x 18	225 x 117 x 88



Test tube racks, anodised aluminium, single platform

03

Thermo
SCIENTIFIC


Catalogue No	Alt. No	Dimensions [w x d x h], mm	Hole diameter, mm	Array
STK-765-010A	RS10	62 x 127 x 28	7.0	5 x 10
STK-765-030R	RS21	207 x 107 x 22	14.2	5 x 10
STK-765-040X	RS31	232 x 98 x 33	17.5	4 x 10
STK-765-050L	RS32	232 x 120 x 33	17.5	5 x 10
STK-765-060Y	RS70	252 x 129 x 29	18.0	5 x 10
STK-765-065V	RS71	252 x 249 x 29	18.0	10 x 10
STK-765-080C	RS80	265 x 133 x 29	19.5	5 x 10
STK-765-100W	RS40	177 x 173 x 29	23.2	6 x 6
STK-765-115J	RS51	208 x 140 x 40	28.2	4 x 6
STK-765-120Q	RS52	208 x 205 x 40	28.2	6 x 6

Metric racks

Catalogue No	Alt. No	Dimensions [w x d x h], mm	Hole diameter, mm	Array
STK-766-010R	MR0-01	187 x 187 x 44	13.2	10 x 10
STK-766-015H	MR002	187 x 280 x 46	13.2	10 x 15
STK-766-020X	MR003	210 x 210 x 61	14.2	10 x 10
STK-766-025E	MR004	210 x 315 x 61	14.2	10 x 15
STK-766-030L	MR005	233 x 350 x 61	17.5	10 x 15
STK-766-035B	MR006	283 x 283 x 61	20.5	10 x 10
STK-766-040Y	MR007	283 x 450 x 61	20.5	10 x 15

Double platform racks

Catalogue No	Alt. No	Dimensions [w x d x h], mm	Hole diameter, mm	Array
STK-765-110T	R155	233 x 120 x 61	17.5	5 x 10



Test tube racks, coldstore aluminium

03

Thermo
SCIENTIFIC

Subdivided into square sections for storage of bottles and tubes. All dimensions nominal.

Catalogue No	No. of sections	Array	Dimensions, mm	Dimensions [w x d x h], mm
STK-760-030B	144	12 x 12	12 x 12 (section)	154 x 154 x 34
STK-760-034Q	288	12 x 12	12 x 12 (section)	307 x 154 x 34
STK-760-070M	81	9 x 9	16 x 16 (section)	154 x 154 x 34
STK-760-078T	198	9 x 22	16 x 16 (section)	384 x 154 x 34
STK-760-080J	64	8 x 8	18 x 18 (section)	154 x 154 x 34
STK-760-090G	36	6 x 6	24 x 24 (section)	154 x 154 x 34
STK-760-098N	90	6 x 15	24 x 24 (section)	384 x 154 x 34
STK-760-130U	25	5 x 5	30 x 30 (section)	154 x 154 x 34
STK-760-150X	50	5 x 10	30 x 30 (section)	307 x 154 x 34
STK-760-170Y	60	5 x 12	30 x 30 (section)	384 x 154 x 34

